



The
University
Of
Sheffield.

NHS
Sheffield Teaching Hospitals
NHS Foundation Trust

yorkshire
cancer
research

POLARIS
Pulmonary, Lung and Respiratory Imaging
Sheffield
Polarised Imaging Systems

MIDL
Lübeck 2021

A hybrid model- and deep learning-based framework for functional lung image synthesis from non-contrast multi-inflation CT


Joshua R. Astley^{1,2}, Alberto M. Biancardi¹, Helen Marshall¹, Guilhem J. Collier¹, Paul J.C. Hughes¹, Jim M. Wild¹, Bilal A. Tahir^{1,2}

¹POLARIS, Department of Infection, Immunity & Cardiovascular Disease and ²Department of Oncology and Metabolism, University of Sheffield, UK


Introduction

- Hyperpolarised gas MRI is a functional lung imaging modality capable of visualising regional ventilation with exquisite detail. However, the modality requires highly specialised equipment and a contrast agent such as Helium-3 (³He) or Xenon-129 (¹²⁹Xe).
- CT ventilation imaging (CTVI) aims to derive images of regional ventilation from multi-inflation CT, acquired during tidal breathing or breath-hold, without exogenous contrast.


E.g. Tahir et al 2018, IJROBP; Tahir et al 2019, Phys Med Biol



³He MRI



¹²⁹Xe MRI



CTVI


Aim

Can deep learning improve upon state-of-the-art model-based CTVI methods?


Methods

Generation of CTVI

1. Acquire inspiratory & expiratory CT
Can be from 4DCT or breath-hold



2. Perform deformable image registration



3. Computational modelling

$$\frac{\Delta V}{V_{exp}} = 1000 \frac{\overline{HU}_{ins} - HU_{exp}}{HU_{exp}(1000 + HU_{ins})}$$


CT^{HU} model



We used four combinations of input channels for the CNN:

- 1) Expiratory CT + inspiratory CT + CTVI
- 2) Expiratory CT + inspiratory CT
- 3) Expiratory CT
- 4) Inspiratory CT

Deep learning methods



Results / Conclusion

Synthetic ventilation generation methods	Spearman's ρ	MSE
	Mean \pm SD	Mean \pm SD
CT ^{HU} model	0.39 \pm 0.18	N/A
DL (expiration CT)	0.41 \pm 0.18	0.032 \pm 0.01
DL (inspiration CT)	0.37 \pm 0.20	0.027 \pm 0.01
DL (expiration CT + inspiration CT)	0.42 \pm 0.18	0.027 \pm 0.01
DL (expiration CT + inspiration CT + CT ^{HU} model)	0.46 \pm 0.16	0.025 \pm 0.01

The hybrid model/DL method achieved a mean Spearman's correlation on 47 scans of 0.46 and significantly outperformed other DL methods as well as CTVI modelling.



Qualitatively, we can see for 3 cases that the hybrid model / DL-based approach is able to accurately replicate defects present in the hyperpolarised gas MRI scans.

